

Isothermal Titration Calorimetry (ITC)


TA Instruments

TAM


Courses

Titration Calorimetry to Measure Binding Interactions

- Universal
- No chemical modification required
- Immobilisation possible, but not necessary
- Equilibrium conditions



Titration calorimetry is the controlled addition (in steps or continuously) of one substance to another, whereby a response in heat flow is measured. Liquid-liquid titration is most common, but solid-liquid (dissolution) and gas-solid (sorption) can also be considered as titrations.



TAM

Courses

The titration experiment


The relation between ΔH and K
(1:1 binding)



$$K = [ML] / [L][M]$$

$$Q = ([ML] / [M]_{\text{total}}) \cdot \Delta H$$

$$Q = (K \cdot [L] \cdot \Delta H) / (1 + K \cdot [L])$$




The Titration ampoule

For studying interaction phenomena;

-in liquid solution


-between solid state samples and gas



Computer controlled
injection pump

TAM

Courses



Reaction
vessel

Typical Application areas

- Interactions of biological macromolecules (e.g. ligand binding, protein-protein interaction)
- Interaction between drugs and release modifiers (e.g. cyclodextrins, PVP)
- Adsorption/binding to particles in suspension
- Dissolution/dilution of liquids
- Gas titration

Binding Characterization

- Binding Affinity, K
- Binding Stoichiometry
- Binding Thermodynamics: ΔG , ΔH , ΔS

$$\Delta G^{\circ} = -RT \ln K_c$$

$$\Delta G^{\circ} = \Delta H^{\circ} - T\Delta S^{\circ}$$

$$\Delta C_p = \frac{d\Delta H^{\circ}}{dT}$$

Titration Calorimetry in Biochemistry/Biophysics

- Emerging as a premier tool for the characterization of biological macromolecules
- General applicability with high accuracy and precision
- Example:
 - Antigen – Antibody
 - Peptide – Protein
 - Lipid – Protein
 - Nucleic Acid – Protein
 - Carbohydrate – Protein
 - Small Molecule/Drug – Protein
 - Protein – Protein
 - Protein – Receptor (soluble and membrane-bound)


TAM

Courses

Rational Drug Design

-Screening of pharmaceutical hits and leads

Hallén, D., Structural Chemistry, Pharmacia, Stockholm





Non-binder: the identical peaks reflects only a dilution process

Binder: the peaks reflects dilution as well as a reaction. The peaks indicate a binding of the drug to a target molecule

Proton Linkage

Hallén, D., Structural Chemistry, Pharmacia & Upjohn, Stockholm

- $M + L + B \rightleftharpoons ML + BH^+$; $\Delta H_{\text{reaction}}$
- $M + L \rightleftharpoons ML + H^+$; $\Delta H_{\text{binding}}$
- $B + H^+ \rightleftharpoons BH^+$; ΔH_{proton}





Binding of IGF-I to the soluble extracellular IGF-I receptor

Binding of a ligand to a macrocyclic compound

Binding of $BaCl_2$ to 18-crown-6-ether

Integrated heat values plotted against ligand/host molecule



The binding process is characterised by the equilibrium constant K_c and the reaction enthalpy ΔH .



$$K=5900 \text{ dm}^3/\text{mol}$$

$$\Delta H = -31.42 \text{ kJ/mol.}$$

TAM

Courses

Binding of a calixarene derivative CATS and alkyl ammonium ions



Thermo-
dynamics

$$\Delta G^{\circ} = -RT \ln K_c$$

$$\Delta G^{\circ} = \Delta H^{\circ} - T\Delta S^{\circ}$$

$$\Delta C_p = \frac{d\Delta H^{\circ}}{dT}$$

Conclusion

The affinity decreases and the reaction becomes more exothermic with increasing chain length

n	K_c (M^{-1})	ΔH (kJ/mol)
3	13090	-16.89
4	10200	-17.94
5	6400	-20.24
6	4000	-20.42
7	2460	-20.86

TAM

Courses

Dr. M. Stödeman, Univ. Lund



IUPAC Test Reaction for Isothermal Titration Calorimetry

- The reaction between 2'CMP and RNaseA has been suggested by IUPAC as a test reaction in titration calorimetry
 - 2'CMP: 2-cytidine monophosphate
 - RNaseA: bovine pancreatic ribonuclease
- 1:1 binding reaction

Dr. Dan Hallén and Dr. Natalia Markova at Biovitrum AB, Stockholm, Sweden.

Experimental procedure and results

- Experimental conditions
 - Conc of RNaseA in vessel: 66.670 $\mu\text{mol/l}$
 - Conc of 2'CMP in syringe: 764.702 $\mu\text{mol/l}$
 - Volume in vessel: 900.000 μl
 - Stirrer rate: 60 rpm
 - 20 mM Potassium Acetate buffer at pH5.5.
- TAM System configuration
 - 1 ml titration ampoule
 - gold stirrer
 - titration reference ampoule
 - 3 μW measuring range
 - Dynamically corrected data



Results

- The derived thermodynamic properties are highly dependent on pH, ionic strength, temperature and RNaseA concentration.

Therefore it is necessary to get all these parameters correct during the sample preparation.

- It was proposed that the standardised test should be performed at high salt concentrations because the enthalpy change is less sensitive to protein concentrations under these conditions.

Ligand binding



$$\frac{ML}{L \cdot M} = K$$

$$L_t = L + ML \quad \Rightarrow \quad ML = L_t - L$$

$$M_t = M + ML \quad \Rightarrow \quad M = M_t - ML = M_t - L_t + L$$

$$q = V \cdot \Delta H \cdot ML$$

$$\frac{1}{V \cdot \Delta H} \frac{dq}{dL_t} = \frac{d}{dL_t} (L_t - L) = 1 - \frac{dL}{dL_t}$$

$$r = K \cdot M_t \quad \Rightarrow \quad M_t = \frac{r}{K}$$


$$X = \frac{L_t}{M_t} \quad \Rightarrow \quad L_t = \frac{X \cdot r}{K}$$

L – ligand (syringe)
M – host molecule (vessel)
K – equilibrium constant
 ΔH – reaction enthalpy

$$\frac{dq}{dL_t} = V \cdot \Delta H \frac{1}{2} \left(1 - \frac{1 + r \cdot (X - 1)}{\sqrt{(r - X \cdot r + 1)^2 + 4X \cdot r}} \right)$$

The left side can be plotted versus $X = L_t/M_t$ for different values of $r = K \cdot M_t$.

Note: the shape of the 'binding curve' depends on the product $K \cdot M_t$. Thus, for a certain K value the concentration of M must be optimized for highest accuracy in the curve fitting procedure.






Overall / Step-wise Reaction Path




General reaction stoichiometry when L binds to M



Overall reaction path



Step-wise reaction path



NOTE: In the Ligand binding analysis program of TAM Assistant the equilibrium constants (β_n) and the reaction enthalpies for the overall reactions (ΔH_i) are determined.

$$\begin{aligned} \beta_1 &= K_1 & \Delta H_{1s} &= \Delta H_1 \\ \beta_2 &= K_1 K_2 & \Delta H_{2s} &= \Delta H_2 - \Delta H_1 \\ \beta_{31} &= K_1 K_2 K_3 & \Delta H_{3s} &= \Delta H_3 - \Delta H_2 \end{aligned}$$

TAM

Courses

Conclusion

- Isothermal Titration Calorimetry (ITC) is a versatile technique that can be used to study a variety of applications
- Some characteristic of the TAM system are;
 - The calorimeter is of heat conduction type (temperature remains constant during the reaction)
 - High dynamic measuring range, i.e up to 3 μW
 - Easy to clean the reaction vessel
 - Different sizes of reaction vessels
 - The reaction vessel can be observed
 - The stirring rate can be varied
 - TAM is flexible allowing the user to developed his own technical solutions

END

TAM

Courses